

3D APPLICATIONS

Cancer Invasion

MATERIAL: TeloCol-6 , TeloCol-4

CELL LINE: MDA-MB-231

CONSTRUCT TYPE: Droplet-in-Droplet


VESSEL TYPE: 96 well-plate

INNER DROPLET: 1 μ L

OUTER DROPLET: 5 μ L

CELLINK PROTOCOL

Understand what causes cancer cells to migrate from diseased to healthy tissue. The droplet-in-droplet method, made possible by 3D biodeposition, enables researchers to create models to understand cancer invasion mechanics. In this proof of concept cancer cells embedded in TeloCol-6 were deposited and then a subsequent droplet of TeloCol-4 was dispensed around to surround the first drop and create a cancerous core.



MDA-MB-231 (mCherry): 4×10^6 cells/mL

TeloCol: 6 mg/mL




FIGURE 1: General workflow for TeloCol-10 biodispersing of cancer cell invasion (droplet-in-droplet) model and analysis capable on systems like the BIO CELLX.




FIGURE 2: Droplet-in-droplet model for cancer cell migration. The live cancer cells in red are embedded in the inner droplet at day 1 and migrated out to the outer layer during the 9 days incubation time.



MORE DETAILS ABOUT THE TUMOR SPHEROIDS
Read more in the full
application note

CONTACT US FOR MORE INFO

www.cellink.com sales@cellink.com | US: +1 (833) CELLINK | EU: +46 31 128 700

CELLINK
A BICO COMPANY