

Multiplexed Evaluation of Soluble Immune-Checkpoint Molecules in Sera from Metastatic NSCLC Patients Receiving Immunotherapy

Imad Tarhoni, MD, PhD, Cristina L. Fhied, MS, David Gerard, MS, and Jeffrey A. Borgia, PhD Departments of Cell & Molecular Medicine, Pathology, and the Rush Biomarker Development Core; Rush University Medical Center, Chicago, IL 60612.

Abstract

Lung cancer remains the primary cause of cancer-related mortalities worldwide, a fact resulting from late diagnosis after widespread dissemination of the disease has occurred. PD-1/PD-L1 directed immune checkpoint inhibition therapy has changed the paradigm of cancer therapy and survivorship in this group with impressive, durable clinical benefits and low toxicity profiles. However, this specific approach is imperfect where only a subset of patients benefit from this therapeutic strategy. Multiple alternative immune checkpoint systems have been described as possible therapeutic options, though their potential impact in lung cancer remains undetermined. In this study, we profiled a panel of 31 immune checkpoint molecules (listed in Table 1), consisting of soluble checkpoint molecules and immune regulators, using the MILLIPLEX® Human Immuno-Oncology Checkpoint Protein Panel 2 with pretreatment sera from 123 advanced-stage Non-Small Cell Lung Cancer patient samples receiving PD-1/PD-L1 targeted therapy. Understanding the dynamics of these molecules and their association with clinical outcomes may help advance research in novel strategies for immune checkpoint inhibition therapy. Results from this evaluation demonstrate that this panel showed excellent analytical performance and identified several soluble checkpoint molecules, including CD40L, APRIL, and CD226, which have a significant (p -values <0.05 for all) association with both progression-free and overall survival rates. Research evaluating soluble immune checkpoints in serum may help identify those likely to benefit from PD-1/PD-L1 immunotherapy and may reveal new immunotherapy targets for future investigation.

Introduction

The immune system is an essential regulator of tumor biology with the potential of supporting or suppressing tumor growth and progression[1]. Harnessing the immune system against cancer is the focus of several therapeutic approaches in cancer research. Immune checkpoint inhibitors (ICIs) are therapeutic monoclonal antibodies designed to disrupt inhibitory signals received by immune checkpoint regulatory molecules and show long-term survival benefits for some patients with metastatic melanoma, Non-Small Cell Lung Cancer (NSCLC), and other tumors[2-4]. However, response to ICIs is not homogenous for all tumors, suggesting heterogeneity of immunoregulatory mechanisms in cancer[5]. The coexistence of multiple checkpoints in the same tumor, for example, imposes a challenge for cancer immunotherapy and could be associated with resistance and treatment failure[6]. Therefore, research to understand these resistance mechanisms may facilitate selecting patients who are most likely to benefit from current therapeutic strategies or indicate new therapeutic targets for future studies. Immune checkpoint molecules are cell-surface proteins expressed on immune cells, mainly T cells, that regulate immune activation by various antigens, including tumor antigens. Immunotherapeutic agents harness the intrinsic immune response against tumor antigens by removing inhibitory mechanisms for T-cell activation upon interaction with antigen-presenting cells[7].

Programmed death receptor-1 (PD-1) and its ligand-1 (PD-L1) are the main targets for current immunotherapy agents approved for NSCLC, the most common subtype of lung cancer. However, the function and potential therapeutic value of numerous other immune checkpoints are not completely understood.

NSCLC is a leading cause of cancer mortality globally and is the second cancer that has benefitted from immune checkpoint therapy, especially in advanced stages[8, 9]. Immunotherapeutic agents that target the PD-1/PD-L1 immune checkpoint pathway have been meteorically incorporated into the standard management of advanced-stage NSCLC after outstanding clinical benefits have been documented in a series of clinical trials. However, overall response remains dismal for these patients and multiple resistance mechanisms have started to emerge, including the expression of alternate immune checkpoint molecules. Here, we evaluate circulating levels of 31 soluble immune checkpoints from patients receiving immunotherapy and correlated the levels to their clinical outcomes as a means to better appreciate the potential involvement of alternate ICI mechanisms in NSCLC.

Table 1: Average intra-assay precision values (expressed as %CV) between replicates. Coefficient of variance (CV) calculated between the replicates of a standard curve, samples, and kit Quality Controls (QCs).

| ANALYTE | MEAN %CV IN STANDARD CURVE | MEAN %CV IN SAMPLES | MEAN %CV IN KIT QCS |
|---------------------|----------------------------|---------------------|---------------------|
| CD40L | 1.9 | 2.3 | 4.0 |
| 4-1BBL/TNFSF9 | 1.6 | 5.8 | 3.2 |
| ARGINASE-1 | 1.9 | 4.3 | 3.1 |
| B7-H2/ICOSL | 2.6 | 2.8 | 1.1 |
| B7-H3/CD276 | 3.4 | 1.9 | 1.4 |
| 5-NT/CD73 | 2.1 | 2.5 | 1.2 |
| B7-H4/VTCN1 | 3.2 | 5.2 | 3.9 |
| APRIL | 3.3 | 2.3 | 0.4 |
| B7-H5/VISTA | 1.8 | 4.5 | 1.6 |
| CD25/IL-2R α | 2.6 | 2.3 | 2.6 |
| B7-H6 | 2.9 | 4.5 | 2.9 |
| CD137/4-1BB | 2.1 | 2.8 | 2.7 |
| GRANZYME B | 2.2 | 3.0 | 1.3 |
| CD226/DNAM-1 | 2.2 | 4.4 | 3.3 |
| CD30/TNFRSF8 | 2.8 | 4.7 | 2.3 |
| E-CADHERIN | 3.0 | 3.2 | 1.6 |
| FGL1/HEPASSOCIN | 1.8 | 1.9 | 1.7 |
| GALACTIN-1 | 2.3 | 1.5 | 1.0 |
| GALACTIN-3 | 2.3 | 1.8 | 2.0 |
| GRANULYCIN | 2.0 | 2.4 | 1.0 |
| IDO1 | 1.8 | 5.0 | 2.8 |
| MICA | 3.6 | 4.4 | 1.0 |
| MICB | 3.4 | 3.8 | 1.7 |
| NECTIN-2 | 2.1 | 1.9 | 1.0 |
| BAFF/BLYS | 2.7 | 1.8 | 1.3 |
| NECTIN-4 | 1.1 | 3.3 | 1.4 |
| OX40/CD134 | 3.0 | 2.5 | 0.7 |
| PVR/CD155 | 4.2 | 2.2 | 4.1 |
| SIGLEC-7 | 2.7 | 4.0 | 2.2 |
| SIGLEC-9 | 2.4 | 3.0 | 1.9 |
| PERFORIN | 3.5 | 1.9 | 1.4 |

Methods

Patient Sample Population:

Serum was collected from 123 patients with advanced-stage NSCLC (IIIb-IV) that failed previous rounds of platinum doublet chemotherapy and before initiation of PD-1/PD-L1 directed immunotherapy (e.g., nivolumab, atezolizumab, pembrolizumab, or durvalumab). Clinical data were collected prospectively after written informed consent was obtained from each subject by the Rush University Cancer Center (RUCC) Biorepository. Clinical outcomes were assigned and calculated from the first immunotherapy administration date as follows: Progression-free survival (PFS) is the time elapsed from the treatment start date until the first tumor progression is encountered. Overall survival (OS) is the time elapsed between the treatment start date and death date or last date of follow-up with the patient known alive. The Rush University Medical Center (RUMC) Institutional Review Board approved all protocols used in this study.

Serum Collection and Storage:

Peripheral blood was collected in 10 mL red-top Vacutainers[®] immediately before initiation of PD-1/PD-L1 directed immunotherapy and processed using standard clinical laboratory methods. A portion of each serum sample used for the protein biomarker evaluations using the MILLIPLEX[®] Luminex[®] based immunoassays was supplemented with 25 μ L/mL of the Mammalian Protease Inhibitor Cocktail (Cat. No. P8340-5ML) and 10 μ L/mL of 0.5 M EDTA to minimize further proteolysis after processing. Aliquots were archived at -80 $^{\circ}$ C until testing, with no specimen tested subjected to greater than two freeze-thaw cycles.

Measurements of Serum Biomarker Levels:


Serum specimens were evaluated with 31 soluble checkpoint molecules and immune regulators using the MILLIPLEX[®] Human Immuno-Oncology Checkpoint Protein Panel 2 (Cat. No. HCKP2-11K). This panel is intended for research use only (RUO), and not for use in diagnostic procedures. This panel consists of the following targets: 4-1BBL/TNFSF9, 5'-NT/CD73 (5'-nucleotidase), APRIL (TNFSF13), Arginase-1, B7-H2/ICOSL (Inducible T cell costimulator ligand), B7-H3/CD276, B7-H4/VTCN1 (V-set domain-containing T-cell activation inhibitor 1), B7-H5/VISTA (V-set immunoregulatory receptor), B7-H6 (Natural killer cell cytotoxicity receptor 3 ligand 1, NCR3LG1), BAFF/BlyS (TNFSF13B), CD25/IL-2R α (Interleukin-2 receptor alpha), CD30/TNFRSF8, CD40L (CD40 ligand), CD137/4-1BB (TNFRSF9), CD226/DNAM-1 (DNAX accessory molecule 1), E-Cadherin, FGL1/Hepassocin (Fibrinogen-like protein 1), Gal-1 (Galectin-1), Gal-3 (Galectin-3), Granulysin, Granzyme B, IDO1 (Indoleamine 2,3-dioxygenase 1), MICA (MHC class

I polypeptide-related sequence A), MICB (MHC class I polypeptide-related sequence B), Nectin-2 (PVRL2, CD112), Nectin-4 (PVRL4), OX40/CD134 (TNFRSF4), Perforin, PVR/CD155 (Poliovirus receptor cell adhesion molecule), Siglec-7 (Sialic acid-binding Ig-like lectin 7, CD328), Siglec-9 (Sialic acid-binding Ig-like lectin 9, CD329). All primary data points were collected on a Luminex® FLEXMAP 3D® system. Analyte concentrations were calculated from a 7-point curve using a five-parametric fit algorithm (xPONENT® v4.0.3 Luminex Corp., Austin, TX). All data met minimum quality control thresholds defined in the kit protocol.

Biomarker Statistical Methods:

Cutoff values in relation to survival parameters were determined using the “survminer” package in R program version 3.4. The association between the baseline biomarkers and the clinical outcomes was calculated using log-rank and Kaplan-Meier analysis using Graphpad v.8.3, and a p-value of 0.05 was used as a cutoff for statistical significance.

Figure 1: Soluble immune checkpoints associated with progression-free and overall survival. Top panel Kaplan-Meier curves showing the association between selected biomarkers and progression-free survival (PFS). The curves represent the cases with baseline markers below or above the cutoff point. **Bottom panel** Kaplan-Meier curves showing the association between selected biomarkers and overall survival (OS). The curves represent the cases with baseline markers below or above the cutoff point.



Results

Analytical Performance:

We quantified 31 analytes associated with immune regulation using the MILLIPLEX® Human Immuno-Oncology Checkpoint Protein Panel 2 in serum specimens representative of the pretreated advanced-stage NSCLC. The assays performed at a wide dynamic range with power of magnitude ranging from three to

four. The precision was averaged at 94% with a high agreement between the samples and the standard curve. Table 1 shows the average coefficient of variation (CV) between the replicates, expressed as %CV (intra-assay precision).

Performance with Patient Cohort:

Based on the cohort we tested, we observed a significant association between baseline soluble immune checkpoint molecules and the clinical outcomes. High baseline serum levels of CD40L, B7-H2/ICOSL, APRIL, B7-H5/VISTA, E-Cadherin, Galectin-1, Galectin-3, IDO1, Nectin-2, OX40/CD134, and Perforin were significantly associated with favorable progression-free survival time p-values <0.05). Kaplan-Meier curves of selected immune checkpoint molecules with cutoffs are depicted in Figure 1 (Top). Hazard ratios and confidence interval statistics 95% CI are shown in detail in Table 2.

Many biomarkers were also associated with overall survival outcome. High baseline levels of CD40L, B7-H2/ICOSL, APRIL, B7-H5/VISTA, E-Cadherin, Galectin-1, Galectin-3, Granulysin, Nectin-2, OX40/CD134, and Perforin were associated with longer overall survival (all p-values <0.05). Lower baseline levels of B7-H3/CD276, 5'-NT/CD73, CD226/DNAM-1, FGL1/Hepassocin, MICA, PVR/CD155, and Siglec-7 were associated with worse overall survival outcome (all p-values <0.05). Hazard ratios and confidence interval statistics 95% CI are shown in detail in Table 3. Figure 1 (Bottom) shows Kaplan-Meier curves of overall survival in relation to baseline cutoff levels in selected immune checkpoint molecules.

Table 2: Soluble immune checkpoints associated with progression-free survival. Association of analytical results with progression-free survival (PFS) with optimized cutoff value indicated, number of cases below ("N≤") and above ("N>"), median PFS below ("LOW") and above ("HIGH") the cutoff, along with the log-rank p-value. Hazard ratios (HR) along with 95% confidence intervals (CI) are also provided.

| TARGET | CUTOFF | N≤ | N> | MEDIAN PFS (LOW/HIGH) | p-VALUE | HR | 95% CI |
|-----------------|-------------|-----|-----|-----------------------|---------|-----|--------------|
| CD40L | 6.2 ng/mL | 60 | 58 | 2.2/6.9 | <0.01 | 0.5 | 0.33 to 0.75 |
| 4-1BBL/TNFSF9 | 62.5 pg/mL | 92 | 36 | 5.4/5.2 | 0.15 | 1.3 | 0.85 to 2.13 |
| ARGINASE-1 | 1.2 ng/mL | 80 | 38 | 5.6/3.4 | 0.17 | 1.3 | 0.85 to 2.09 |
| B7-H2/ICOSL | 7.3 ng/mL | 24 | 94 | 4.1/9.3 | 0.04 | 0.5 | 0.35 to 0.91 |
| B7-H3/CD276 | 991.0 pg/mL | 11 | 107 | 16.3/4.1 | 0.04 | 2.2 | 1.25 to 4.11 |
| 5'-NT/CD73 | 1.9 ng/mL | 27 | 91 | 14.3/3.9 | 0.03 | 1.7 | 1.11 to 2.70 |
| B7-H4/VTCN1 | 484.7 pg/mL | 58 | 65 | 12.6/3.6 | 0.07 | 1.5 | 1.00 to 2.53 |
| APRIL | 2.6 ng/mL | 46 | 72 | 2.2/7.0 | 0.01 | 0.6 | 0.40 to 0.95 |
| B7-H5/VISTA | 90.6 pg/mL | 105 | 13 | 3.9/18.2 | 0.02 | 0.4 | 0.26 to 0.77 |
| CD25/IL-2Rα | 1.6 ng/mL | 17 | 101 | 2.1/5.6 | 0.14 | 0.6 | 0.35 to 1.26 |
| B7-H6 | 82.8 pg/mL | 102 | 16 | 4.1/9.2 | 0.05 | 0.5 | 0.32 to 0.91 |
| CD137/4-1BB | 15.3 pg/mL | 54 | 64 | 7.0/4.0 | 0.13 | 1.3 | 0.90 to 2.03 |
| GRANZYME B | 3.6 pg/mL | 17 | 101 | 10.2/4.8 | 0.06 | 1.8 | 1.07 to 3.09 |
| CD226/DNAM-1 | 95.3 pg/mL | 33 | 85 | 8.6/3.3 | 0.01 | 1.7 | 1.16 to 2.72 |
| CD30/TNFRSF8 | 13.8 pg/mL | 36 | 82 | 6.0/3.9 | 0.12 | 1.4 | 0.92 to 2.15 |
| E-CADHERIN | 25.7 ng/mL | 14 | 104 | 1.5/5.6 | 0.02 | 0.5 | 0.24 to 1.11 |
| FGL1/HEPASSOCIN | 90.7 ng/mL | 83 | 39 | 7.5/2.3 | <0.01 | 1.9 | 1.19 to 3.21 |
| GALECTIN-1 | 22.9 ng/mL | 75 | 43 | 3.7/7.5 | 0.02 | 0.6 | 0.40 to 0.90 |
| GALECTIN-3 | 7.2 ng/mL | 73 | 45 | 3.7/10.0 | 0.01 | 0.5 | 0.38 to 0.86 |
| GRANULYSIN | 786.3 pg/mL | 28 | 90 | 2.1/6.1 | 0.01 | 0.5 | 0.34 to 0.99 |
| IDO1 | 7.7 pg/mL | 12 | 106 | 2.1/5.6 | 0.03 | 0.5 | 0.23 to 1.19 |
| MICA | 7.7 pg/mL | 11 | 107 | 28.6/4.1 | 0.03 | 2.3 | 1.28 to 4.17 |
| MICB | 30.7 pg/mL | 11 | 107 | 22.9/4.1 | 0.05 | 2.1 | 1.17 to 3.93 |
| NECTIN-2 | 1.4 ng/mL | 62 | 56 | 3.9/6.2 | 0.01 | 0.5 | 0.39 to 0.89 |
| BAFF/BLYS | 555.9 pg/mL | 11 | 107 | 16.3/5.2 | 0.24 | 1.5 | 0.82 to 2.97 |
| NECTIN-4 | 307.4 pg/mL | 101 | 17 | 4.1/9.3 | 0.11 | 0.6 | 0.35 to 1.02 |
| OX40/CD134 | 142.9 pg/mL | 54 | 64 | 2.9/6.2 | 0.04 | 0.6 | 0.44 to 1.00 |
| PVR/CD155 | 105.7 pg/mL | 41 | 77 | 9.2/3.3 | 0.03 | 1.5 | 1.04 to 2.38 |
| SIGLEC-7 | 5.2 pg/mL | 43 | 75 | 10.2/4.1 | 0.02 | 1.5 | 1.06 to 2.40 |
| SIGLEC-9 | 4.5 pg/mL | 17 | 101 | 4.1/5.2 | 0.26 | 1.4 | 0.82 to 2.47 |
| PERFORIN | 4.9 ng/mL | 36 | 82 | 2.0/7.0 | <0.01 | 0.4 | 0.29 to 0.80 |

Table 3: Soluble immune checkpoints associated with overall survival. Association of analytical results with overall survival (OS) with optimized cutoff value indicated, number of cases below (“N≤”) and above (“N>”), median PFS below (“LOW”) and above (“HIGH”) the cutoff, along with the log-rank p-value. Hazard ratios (HR) along with 95% confidence intervals (CI) are also provided.

| TARGET | CUTOFF | N≤ | N> | MEDIAN PFS (LOW/HIGH) | p-VALUE | HR | 95% CI |
|-----------------|-------------|-----|-----|-----------------------|---------|-----|--------------|
| CD40L | 6.2 ng/mL | 63 | 60 | 6.6/23.7 | <0.01 | 0.5 | 0.34 to 0.83 |
| 4-1BBL/TNFSF9 | 62.5 pg/mL | 95 | 38 | 18.1/7.6 | 0.14 | 1.3 | 0.85 to 2.24 |
| ARGINASE-1 | 864.5 pg/mL | 66 | 57 | 24.2/10 | 0.02 | 1.7 | 1.08 to 2.58 |
| B7-H2/ICOSL | 5.8 ng/mL | 91 | 32 | 11.4/45.0 | <0.01 | 0.4 | 0.28 to 0.71 |
| B7-H3/CD276 | 995.0 pg/mL | 13 | 110 | 25.7/13.5 | 0.25 | 1.6 | 0.81 to 2.98 |
| 5'-NT/CD73 | 2.6 ng/mL | 44 | 79 | 22.3/13.4 | 0.17 | 1.4 | 0.88 to 2.13 |
| B7-H4/VTCN1 | 164.0 pg/mL | 29 | 94 | 24.7/10.1 | 0.07 | 2.0 | 1.24 to 3.15 |
| APRIL | 2.4 ng/mL | 46 | 77 | 6.6/22.3 | 0.03 | 0.6 | 0.38 to 0.99 |
| B7-H5/VISTA | 90.6 pg/mL | 106 | 17 | 13.4/N/A | 0.05 | 0.5 | 0.26 to 0.84 |
| CD25/IL-2Ra | 411.6 pg/mL | 17 | 106 | 4.3/18.0 | 0.04 | 0.6 | 0.27 to 1.14 |
| B7-H6 | 103.3 pg/mL | 107 | 16 | 13.5/N/A | 0.05 | 0.5 | 0.32 to 0.91 |
| CD137/4-1BB | 10.0 pg/mL | 26 | 97 | 8.7/18.2 | 0.09 | 0.6 | 0.35 to 1.18 |
| GRANZYME B | 3.6 pg/mL | 80 | 43 | 13.0/23.4 | 0.11 | 0.6 | 0.44 to 1.07 |
| CD226/DNAM-1 | 1.2 ng/mL | 55 | 68 | 24.0/10.7 | 0.04 | 1.5 | 1.02 to 2.42 |
| CD30/TNFRSF8 | 28.3 pg/mL | 79 | 44 | 18.2/9.8 | 0.13 | 1.4 | 0.88 to 2.21 |
| E-CADHERIN | 25.7 ng/mL | 14 | 104 | 1.5/5.6 | <0.01 | 0.5 | 0.15 to 0.98 |
| FGL1/HEPASSOCIN | 90.7 ng/mL | 83 | 39 | 7.5/2.3 | <0.01 | 1.9 | 1.19 to 3.21 |
| GALECTIN-1 | 25.2 ng/mL | 98 | 25 | 3.7/7.5 | 0.02 | 0.6 | 0.40 to 0.90 |
| GALECTIN-3 | 9.5 ng/mL | 78 | 45 | 10.7/23.6 | 0.02 | 0.6 | 0.38 to 0.91 |
| GRANULYSIN | 786.3 pg/mL | 29 | 94 | 6.7/21.0 | 0.02 | 0.6 | 0.32 to 0.98 |
| IDO1 | 17.3 pg/mL | 14 | 109 | 11.4/16.1 | 0.06 | 0.6 | 0.25 to 1.24 |
| MICA | 9.4 pg/mL | 16 | 107 | 14.1/14.4 | 0.55 | 1.2 | 0.64 to 2.33 |
| MICB | 228.5 pg/mL | 75 | 48 | 11.4/21.0 | 0.05 | 0.6 | 0.41 to 0.99 |
| NECTIN-2 | 1.2 ng/mL | 33 | 90 | 8.7/19.6 | 0.02 | 0.6 | 0.35 to 0.99 |
| BAFF/BLYS | 628.3 pg/mL | 19 | 104 | 22.0/14.2 | 0.17 | 1.6 | 0.89 to 2.87 |
| NECTIN-4 | 307.4 pg/mL | 61 | 62 | 9.5/18.1 | 0.09 | 0.7 | 0.44 to 1.06 |
| OX40/CD134 | 142.9 pg/mL | 56 | 67 | 9.0/22.0 | 0.02 | 0.6 | 0.39 to 0.94 |
| PVR/CD155 | 93.2 ng/mL | 33 | 90 | 25.7/13.0 | 0.05 | 1.7 | 1.04 to 2.62 |
| SIGLEC-7 | 5.3 pg/mL | 74 | 49 | 18.2/10.1 | 0.19 | 1.3 | 0.85 to 2.08 |
| SIGLEC-9 | 4.5 pg/mL | 17 | 106 | 37.9/14.0 | 0.17 | 1.6 | 0.90 to 3.00 |
| PERFORIN | 4.9 ng/mL | 34 | 89 | 4.0/22.0 | <0.01 | 0.5 | 0.27 to 0.82 |

Conclusions

Soluble immune checkpoint molecules can be evaluated using the MILLIPLEX® Human Immuno-Oncology checkpoint Protein Panel 2 (Cat. No. HCKP2-11K). Moreover, we showed a significant association between the circulating immune checkpoints molecules and the response pattern in NSCLC samples receiving anti-PD-1/PD-L1 therapy.

Discover more at SigmaAldrich.com/milliplex

References

- Schreiber, R.D., L.J. Old, and M.J. Smyth, Cancer immunoediting: integrating immunity's roles in cancer suppression and promotion. *Science*, 2011. 331(6024): p. 1565-70.
- Dafni, U., et al., Immune checkpoint inhibitors, alone or in combination with chemotherapy, as first-line treatment for advanced non-small cell lung cancer. A systematic review and network meta-analysis. *Lung Cancer*, 2019. 134: p. 127-140.
- Keenan, B.P., L. Fong, and R.K. Kelley, Immunotherapy in hepatocellular carcinoma: the complex interface between inflammation, fibrosis, and the immune response. *J Immunother Cancer*, 2019. 7(1): p. 267.
- Shen, X. and B. Zhao, Efficacy of PD-1 or PD-L1 inhibitors and PD-L1 expression status in cancer: meta-analysis. *BMJ*, 2018. 362: p. k3529.
- Li, X., et al., Lessons learned from the blockade of immune checkpoints in cancer immunotherapy. *J Hematol Oncol*, 2018. 11(1): p. 31.
- Li, J., L. Ni, and C. Dong, Immune checkpoint receptors in cancer: redundant by design? *Curr Opin Immunol*, 2017. 45: p. 37-42.
- Wei, S.C., C.R. Duffy, and J.P. Allison, Fundamental Mechanisms of Immune Checkpoint Blockade Therapy. *Cancer Discov*, 2018. 8(9): p. 1069-1086.
- Hellmann, M.D., et al., Nivolumab plus Ipilimumab in Advanced Non-Small-Cell Lung Cancer. *N Engl J Med*, 2019.
- Siegel, R.L., K.D. Miller, and A. Jemal, Cancer statistics, 2020. *CA Cancer J Clin*, 2020. 70(1): p. 7-30.

For Research Use Only. Not For Use In Diagnostic Procedures.

Vacutainers® is a trademark of Becton, Dickinson and Company.

To place an order or receive technical assistance

In Europe, please call Customer Service:

| | |
|-----------------------|-------------------------------|
| France: 0825 045 645 | Spain: 901 516 645 Option 1 |
| Germany: 069 86798021 | Switzerland: 0848 645 645 |
| Italy: 848 845 645 | United Kingdom: 0870 900 4645 |

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: [MerckMillipore.com/offices](https://www.MerckMillipore.com/offices)

For Technical Service visit: [MerckMillipore.com/techservice](https://www.MerckMillipore.com/techservice)

[MerckMillipore.com](https://www.MerckMillipore.com)

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

